

## Product datasheet for **HK300890**

### **MIR3913-2 Human qPCR Template Standard (MI0016418)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>Forward Sequence:</b>	TTTGGGACTGATCTTGATGTC
<b>Reverse Sequence:</b>	GAACATGTCTGCGTATCTC
<b>ACCN:</b>	MIMAT0018187
<b>Synonyms:</b>	hsa-miR-3913-5p; MI0016418; MIMAT0018187; MIR3913-2 hsa-mir-3913-2
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
<b>Quality Control:</b>	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
<b>Storage:</b>	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)