

Product datasheet for **HK300927**

MIR4259 Human qPCR Template Standard (MI0015858)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	CAGTTGGGTCTAGGGTCAG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0016880
Synonyms:	hsa-miR-4259; MI0015858; MIMAT0016880; MIR4259 hsa-mir-4259
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.

