

Product datasheet for **HK300931**

MIR4269 Human qPCR Template Standard (MI0015875)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: GCAGGCACAGACAGCCCTG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0016897

Synonyms: hsa-miR-4269; MI0015875; MIMAT0016897; MIR4269 hsa-mir-4269

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)