

Product datasheet for **HK300941**

MIR4322 Human qPCR Template Standard (MI0015851)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: CTGTGGGCTCAGCGCGTG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0016873

Synonyms: hsa-miR-4322; MI0015851; MIMAT0016873; MIR4322 hsa-mir-4322

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)