

Product datasheet for **HK300963**

MIR4445 Human qPCR Template Standard (MI0016788)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: AGATTGTTTCTTTGCCGTG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0018963

Synonyms: hsa-miR-4445-5p; MI0016788; MIMAT0018963; MIR4445 hsa-mir-4445

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)