

## Product datasheet for **HK301013**

### **MIR4521 Human qPCR Template Standard (MI0016887)**

#### **Product data:**

**Product Type:** qPCR Template Standards

**Forward Sequence:** GCTAAGGAAGTCCTGTGCTC

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0019058

**Synonyms:** hsa-miR-4521; MI0016887; MIMAT0019058; MIR4521 hsa-mir-4521

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)