

Product datasheet for **HK301016**

MIR4524a Human qPCR Template Standard (MI0016891)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: ATAGCAGCATGAACCTGTC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0019062

Synonyms: hsa-miR-4524a-5p; MI0016891; MIMAT0019062; MIR4524a hsa-mir-4524a

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)