

## Product datasheet for **HK301026**

### **MIR4538 Human qPCR Template Standard (MI0016909)**

#### **Product data:**

**Product Type:** qPCR Template Standards

**Forward Sequence:** GAGCTTGGATGAGCTGGGC

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0019081

**Synonyms:** hsa-miR-4538; MI0016909; MIMAT0019081; MIR4538 hsa-mir-4538

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)