

Product datasheet for **HK301033**

MIR4636 Human qPCR Template Standard (MI0017263)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	AACTCGTGTTCAAAGCCTTT
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0019693
Synonyms:	hsa-miR-4636; MI0017263; MIMAT0019693; MIR4636 hsa-mir-4636
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)