

Product datasheet for **HK301037**

MIR4640 Human qPCR Template Standard (MI0017267)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TGGGCCAGGGAGCAGCTGGT

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0019699

Synonyms: hsa-miR-4640-5p; MI0017267; MIMAT0019699; MIR4640 hsa-mir-4640

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)