

## Product datasheet for **HK301040**

### **MIR4644 Human qPCR Template Standard (MI0017271)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>Forward Sequence:</b>	TGGAGAGAGAAAAGAGACAGA
<b>Reverse Sequence:</b>	GAACATGTCTGCGTATCTC
<b>ACCN:</b>	MIMAT0019704
<b>Synonyms:</b>	hsa-miR-4644; MI0017271; MIMAT0019704; MIR4644 hsa-mir-4644
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
<b>Quality Control:</b>	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
<b>Storage:</b>	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)