

## Product datasheet for **HK301045**

### **MIR4651 Human qPCR Template Standard (MI0017279)**

#### **Product data:**

**Product Type:** qPCR Template Standards

**Forward Sequence:** CGGGGTGGGTGAGGTCGG

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0019715

**Synonyms:** hsa-miR-4651; MI0017279; MIMAT0019715; MIR4651 hsa-mir-4651

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)