

Product datasheet for **HK301072**

MIR4680 Human qPCR Template Standard (MI0017312)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	AGAACTCTTGCACTCTTAGATG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0019764
Synonyms:	hsa-miR-4680-5p; MI0017312; MIMAT0019764; MIR4680 hsa-mir-4680
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.

