

Product datasheet for **HK301096**

MIR4706 Human qPCR Template Standard (MI0017339)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: AGCGGGGAGGAAGTGGGCGCT

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0019806

Synonyms: hsa-miR-4706; MI0017339; MIMAT0019806; MIR4706 hsa-mir-4706

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)