

Product datasheet for **HK301097**

MIR4707 Human qPCR Template Standard (MI0017340)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	GCCCCGGCGCGGGCGGGTTC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0019807
Synonyms:	hsa-miR-4707-5p; MI0017340; MIMAT0019807; MIR4707 hsa-mir-4707
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)