

Product datasheet for **HK301102**

MIR4713 Human qPCR Template Standard (MI0017347)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TTCTCCCACTACCAGGCTC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0019820
Synonyms:	hsa-miR-4713-5p; MI0017347; MIMAT0019820; MIR4713 hsa-mir-4713
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)