

## Product datasheet for **HK301112**

### MIR4725 Human qPCR Template Standard (MI0017362)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	AGACCCTGCAGCCTTCCC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0019843
Synonyms:	hsa-miR-4725-5p; MI0017362; MIMAT0019843; MIR4725 hsa-mir-4725
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)