

## Product datasheet for **HK301132**

### **MIR4746 Human qPCR Template Standard (MI0017385)**

#### **Product data:**

|                          |                                                                                                                                                                                                       |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b>     | qPCR Template Standards                                                                                                                                                                               |
| <b>Forward Sequence:</b> | CCGGTCCCAGGAGAACCTG                                                                                                                                                                                   |
| <b>Reverse Sequence:</b> | GAACATGTCTGCGTATCTC                                                                                                                                                                                   |
| <b>ACCN:</b>             | MIMAT0019880                                                                                                                                                                                          |
| <b>Synonyms:</b>         | hsa-miR-4746-5p; MI0017385; MIMAT0019880; MIR4746 hsa-mir-4746                                                                                                                                        |
| <b>Components:</b>       | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions) |
| <b>Quality Control:</b>  | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment                                                                                                     |
| <b>Storage:</b>          | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.                                                                                                   |

