

Product datasheet for **HK301140**

MIR4754 Human qPCR Template Standard (MI0017394)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: ATGCGGACCTGGGTTAGCG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0019894

Synonyms: hsa-miR-4754; MI0017394; MIMAT0019894; MIR4754 hsa-mir-4754

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)