

Product datasheet for **HK301176**

MIR4999 Human qPCR Template Standard (MI0017865)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TGCTGTATTGTCAGGTAGTG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0021017
Synonyms:	hsa-miR-4999-5p; MI0017865; MIMAT0021017; MIR4999 hsa-mir-4999
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.

