

Product datasheet for **HK301210**

MIR5197 Human qPCR Template Standard (MI0018176)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Template Standards |
| Forward Sequence: | CAATGGCACAAACTCATTCTTG |
| Reverse Sequence: | GAACATGTCTGCGTATCTC |
| ACCN: | MIMAT0021130 |
| Synonyms: | hsa-miR-5197-5p; MI0018176; MIMAT0021130; MIR5197 hsa-mir-5197 |
| Components: | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)