

Product datasheet for **HK301212**

MIR548aa-1 Human qPCR Template Standard (MI0016689)

Product data:

| | |
|--------------------------|---|
| Product Type: | qPCR Template Standards |
| Forward Sequence: | AAAAACCAATTACTTTTGAC |
| Reverse Sequence: | GAACATGTCTGCGTATCTC |
| ACCN: | MIMAT0018447 |
| Synonyms: | hsa-miR-548aa; MI0016689; MIMAT0018447; MIR548aa-1 hsa-mir-548aa-1 |
| Components: | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)