

Product datasheet for **HK301223**

MIR548w Human qPCR Template Standard (MI0014222)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: AAAAGTAACTGCGGTTTTTGCC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0015060

Synonyms: hsa-miR-548w; MI0014222; MIMAT0015060; MIR548w hsa-mir-548w

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)