

Product datasheet for **HK301227**

MIR550b-2 Human qPCR Template Standard (MI0016687)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: ATGTGCCTGAGGGAGTAAGA

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0022737

Synonyms: hsa-miR-550b-2-5p; MI0016687; MIMAT0022737; MIR550b-2 hsa-mir-550b-2

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)