

## Product datasheet for **HK301276**

### MIR6087 Human qPCR Template Standard (MI0020364)

#### Product data:

|                   |                                                                                                                                                                                                       |
|-------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Template Standards                                                                                                                                                                               |
| Forward Sequence: | TGAGGCGGGGGGCGAG                                                                                                                                                                                      |
| Reverse Sequence: | GAACATGTCTGCGTATCTC                                                                                                                                                                                   |
| ACCN:             | MIMAT0023712                                                                                                                                                                                          |
| Synonyms:         | hsa-miR-6087; MI0020364; MIMAT0023712; MIR6087 hsa-mir-6087                                                                                                                                           |
| Components:       | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions) |
| Quality Control:  | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment                                                                                                     |
| Storage:          | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.                                                                                                   |

