

Product datasheet for **HK301277**

MIR6088 Human qPCR Template Standard (MI0020365)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: AGAGATGAAGCGGGGGG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0023713

Synonyms: hsa-miR-6088; MI0020365; MIMAT0023713; MIR6088 hsa-mir-6088

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)