

## Product datasheet for **HK301278**

### MIR6089-1 Human qPCR Template Standard (MI0020366)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	GGAGGCCGGGTGGGGCGGGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0023714
Synonyms:	hsa-miR-6089; MI0020366; MIMAT0023714; MIR6089-1 hsa-mir-6089-1
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)