

Product datasheet for **HK301286**

MIR6133 Human qPCR Template Standard (MI0021278)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: TGAGGGAGGAGGTTGGG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0024617

Synonyms: hsa-miR-6133; MI0021278; MIMAT0024617; MIR6133 hsa-mir-6133

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)