

## Product datasheet for **HK301301**

### MIR6511a-2 Human qPCR Template Standard (MI0023564)

#### Product data:

**Product Type:** qPCR Template Standards

**Forward Sequence:** CAGGCAGAAGTGGGGCTGA

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0025478

**Synonyms:** hsa-miR-6511a-5p; MI0023564; MIMAT0025478; MIR6511a-2 hsa-mir-6511a-2

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)