

Product datasheet for **HK301313**

MIR6715b Human qPCR Template Standard (MI0022549)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: ACAGGCACGACTGGTTTGG

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0025842

Synonyms: hsa-miR-6715b-5p; MI0022549; MIMAT0025842; MIR6715b hsa-mir-6715b

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)