

Product datasheet for **HK301316**

MIR6718 Human qPCR Template Standard (MI0022553)

Product data:

| | |
|--------------------------|---|
| Product Type: | qPCR Template Standards |
| Forward Sequence: | TAGTGGTCAGAGGGCTTAT |
| Reverse Sequence: | GAACATGTCTGCGTATCTC |
| ACCN: | MIMAT0025849 |
| Synonyms: | hsa-miR-6718-5p; MI0022553; MIMAT0025849; MIR6718 hsa-mir-6718 |
| Components: | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)