

## Product datasheet for **HK301320**

### MIR6722 Human qPCR Template Standard (MI0022557)

#### Product data:

**Product Type:** qPCR Template Standards

**Forward Sequence:** AGGCGCACCCGACCACAT

**Reverse Sequence:** GAACATGTCTGCGTATCTC

**ACCN:** MIMAT0025853

**Synonyms:** hsa-miR-6722-5p; MI0022557; MIMAT0025853; MIR6722 hsa-mir-6722

**Components:** 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10<sup>7</sup> copies (100 reactions)

**Quality Control:** Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

**Storage:** The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)