

## Product datasheet for **HK301354**

### MIR6758 Human qPCR Template Standard (MI0022603)

#### Product data:

|                   |   |
|-------------------|---|
| Product Type:     | qPCR Template Standards   |
| Forward Sequence: | TAGAGAGGGGAAGGATGTGA  |
| Reverse Sequence: | GAACATGTCTGCGTATCTC   |
| ACCN:             | MIMAT0027416  |
| Synonyms:         | hsa-miR-6758-5p; MI0022603; MIMAT0027416; MIR6758 hsa-mir-6758  |
| Components:       | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions) |
| Quality Control:  | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment   |
| Storage:          | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.   |



[View online »](#)