

Product datasheet for **HK301359**

MIR6762 Human qPCR Template Standard (MI0022607)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: CGGGGCCATGGAGCAGCCT

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0027424

Synonyms: hsa-miR-6762-5p; MI0022607; MIMAT0027424; MIR6762 hsa-mir-6762

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)