

Product datasheet for **HK301361**

MIR6764 Human qPCR Template Standard (MI0022609)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TCCCAGGGTCTGGTCAGAG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0027428
Synonyms:	hsa-miR-6764-5p; MI0022609; MIMAT0027428; MIR6764 hsa-mir-6764
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.

