

Product datasheet for **HK301372**

MIR6772 Human qPCR Template Standard (MI0022617)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TGGGTGTAGGCTGGAGCTG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0027444
Synonyms:	hsa-miR-6772-5p; MI0022617; MIMAT0027444; MIR6772 hsa-mir-6772
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)