

Product datasheet for **HK301377**

MIR6777 Human qPCR Template Standard (MI0022622)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	ACGGGGAGTCAGGCAGTG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0027454
Synonyms:	hsa-miR-6777-5p; MI0022622; MIMAT0027454; MIR6777 hsa-mir-6777
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)