

Product datasheet for **HK301381**

MIR6780b Human qPCR Template Standard (MI0022681)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TGGGGAAGGCTTGGCAGGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0027572
Synonyms:	hsa-miR-6780b-5p; MI0022681; MIMAT0027572; MIR6780b hsa-mir-6780b
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.

