

Product datasheet for **HK301412**

MIR6811 Human qPCR Template Standard (MI0022656)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	ATGCAGGCCTGTGTACAGC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0027522
Synonyms:	hsa-miR-6811-5p; MI0022656; MIMAT0027522; MIR6811 hsa-mir-6811
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)