

## Product datasheet for **HK301418**

### MIR6817 Human qPCR Template Standard (MI0022662)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TCTGCCATAGGAAGCTTGGAG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0027534
Synonyms:	hsa-miR-6817-5p; MI0022662; MIMAT0027534; MIR6817 hsa-mir-6817
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)