

## Product datasheet for **HK301427**

### **MIR6826 Human qPCR Template Standard (MI0022671)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>Forward Sequence:</b>	TCAATAGGAAAGAGGTGGGAC
<b>Reverse Sequence:</b>	GAACATGTCTGCGTATCTC
<b>ACCN:</b>	MIMAT0027552
<b>Synonyms:</b>	hsa-miR-6826-5p; MI0022671; MIMAT0027552; MIR6826 hsa-mir-6826
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
<b>Quality Control:</b>	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
<b>Storage:</b>	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)