

Product datasheet for **HK301434**

MIR6833 Human qPCR Template Standard (MI0022678)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: GTGTGGAAGATGGGAGGAGA

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0027566

Synonyms: hsa-miR-6833-5p; MI0022678; MIMAT0027566; MIR6833 hsa-mir-6833

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)