

## Product datasheet for **HK301443**

### MIR6842 Human qPCR Template Standard (MI0022688)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TGGGGTGGTCTCTAGCC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0027586
Synonyms:	hsa-miR-6842-5p; MI0022688; MIMAT0027586; MIR6842 hsa-mir-6842
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)