

Product datasheet for **HK301463**

MIR6861 Human qPCR Template Standard (MI0022708)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: ACTGGGTAGGTGGGGCTC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0027623

Synonyms: hsa-miR-6861-5p; MI0022708; MIMAT0027623; MIR6861 hsa-mir-6861

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)