

## Product datasheet for **HK301464**

### MIR6862-1 Human qPCR Template Standard (MI0022709)

#### Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	CGGGCATGCTGGGAGAGAC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0027625
Synonyms:	hsa-miR-6862-5p; MI0022709; MIMAT0027625; MIR6862-1 hsa-mir-6862-1
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 <sup>7</sup> copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)