

Product datasheet for **HK301497**

MIR6894 Human qPCR Template Standard (MI0022741)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	AGGAGGATGGAGAGCTGGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0027688
Synonyms:	hsa-miR-6894-5p; MI0022741; MIMAT0027688; MIR6894 hsa-mir-6894
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.

