

Product datasheet for **HK301498**

MIR6895 Human qPCR Template Standard (MI0022742)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Template Standards |
| Forward Sequence: | CAGGGCCAGGCACAGAGTA |
| Reverse Sequence: | GAACATGTCTGCGTATCTC |
| ACCN: | MIMAT0027690 |
| Synonyms: | hsa-miR-6895-5p; MI0022742; MIMAT0027690; MIR6895 hsa-mir-6895 |
| Components: | 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions) |
| Quality Control: | Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment |
| Storage: | The primer mix and template standard are stable for one year from date of shipping. Store at -20°C. |



[View online »](#)