

Product datasheet for **HK301500**

MIR7107 Human qPCR Template Standard (MI0022958)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	TCGGCCTGGGGAGGAGGAA
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0028111
Synonyms:	hsa-miR-7107-5p; MI0022958; MIMAT0028111; MIR7107 hsa-mir-7107
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)