

Product datasheet for **HK301524**

MIR7515 Human qPCR Template Standard (MI0024354)

Product data:

Product Type: qPCR Template Standards

Forward Sequence: AGAAGGGAAGATGGTGAC

Reverse Sequence: GAACATGTCTGCGTATCTC

ACCN: MIMAT0029310

Synonyms: hsa-miR-7515; MI0024354; MIMAT0029310; MIR7515 hsa-mir-7515

Components: 1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10⁷ copies (100 reactions)

Quality Control: Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment

Storage: The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.



[View online »](#)