

Product datasheet for **HK301542**

MIR7849 Human qPCR Template Standard (MI0025519)

Product data:

Product Type:	qPCR Template Standards
Forward Sequence:	GACAATTGTTGATCTTGGGC
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0030424
Synonyms:	hsa-miR-7849-3p; MI0025519; MIMAT0030424; MIR7849 hsa-mir-7849
Components:	1 vial of lyophilized SybGREEN qPCR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions), 1 vial lyophilized qPCR miRNA template standard: 50 X 10 ⁷ copies (100 reactions)
Quality Control:	Each standard is generated using qSTAR miRNA primer pairs and sequence-verified prior to shipment
Storage:	The primer mix and template standard are stable for one year from date of shipping. Store at -20°C.

